



- XmnI(124)
- BstBI(128)
- SpeI(132)
- KpnI(142)
- BamHI(144)
- EcoRI(150)
- NruI(156)
- AsiSI(160)
- NotI(163)
- XhoI(171)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	N-His

B01cl MCS region

XmnI
BstBI
SpeI
KpnI
BamHI

...ATGTCGTACTACCATCACCATCACCATCACGAAGGAGTTTCGAACTAGTGGTACCGGATCC
 M S Y Y H H H H H H E G V R T S G T G S

EcoRI
AsiSI
NotI
XhoI

GAATTCGCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 E F A I A A A L E H H H H...