



- XmnI(781)
- BstBI(785)
- SpeI(789)
- KpnI(799)
- BamHI(801)
- EcoRI(807)
- NruI(813)
- AsiSI(817)
- NotI(820)
- XhoI(828)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	GST

B04cl MCS region

XmnI
BstBI
SpeI
KpnI
BamHI
EcoRI

...GTTCCGCGTGATGACGATGACAAGGAAGGAGTTCGA^{ACTAGTGGTACCGGATCCGAATTC}
 GST...V P R **D D D D K** E G V R T S G T G S E F
Enterokinase site ↗

AsiSI
NotI
XhoI

GCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 A I A A A L E H H H H ...