



- XmnI(1379)
- BstBI(1383)
- SpeI(1387)
- KpnI(1397)
- BamHI(1399)
- EcoRI(1405)
- NruI(1411)
- AsiSI(1415)
- NotI(1418)
- XhoI(1426)
- HindIII(1432)

Vector Features

Promoter	Tac
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	MBP

B07cl MCS region

XmnI BstBI SpeI KpnI BamHI
 ...GGGTCGGGCGAAAACCTGTATTTTCAGGGCGAAGGAGTTTCGAACTAGTGGTACCGGATCC
 MBP...G **S G E N L Y F Q G** E G V R T S G T G S
Tev Protease site ↑

EcoRI AsiSI NotI XhoI HindIII
 GAATTCGCGATCGCGGCCGCACTCGAGAAGCTTGGC...
 E F A I A A A L E K L G...