



- XmnI(826)
- BstBI(830)
- KpnI(844)
- BamHI(846)
- EcoRI(852)
- NruI(848)
- AsiSI(862)
- NotI(865)
- XhoI(873)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	His GST

B20cl MCS region

XmnI
BstBI
SpeI
KpnI
BamHI

...CGTTCGGGCGAAAACCTGTATTTTCAGGGCGAAGGAGTTCGAACTAGTGGTACCGGATCC
 GST...R **S G E N L Y F Q G** E G V R T S G T G S
Tev Protease site ↑

EcoRI
AsiSI
NotI
XhoI

GAATTCGCGATCGCGGCCGCACTCGAGCACCACCAC...
 E F A I A A A L E H H H H...