



- XmnI(469)
- BstBI(473)
- SpeI(477)
- KpnI(487)
- BamHI(489)
- EcoRI(495)
- NruI(501)
- AsiSI(505)
- NotI(508)
- XhoI(516)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	HisSUMOAvi

B40cl MCS region

XmnI BstBI SpeI KpnI BamHI EcoRI
 ...GCTCAGAAAATCGAATGGCACGAAGAAAAAGTTTCGAACTAGTGGTACCGGATCCGAATTC
 Avi ...A Q K I E W H E E K V R T S G T G S E F

AsiSI NotI XhoI
 GCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 A I A A A L E H H H H...