



- XmnI(477)
- BstBI(481)
- SpeI(485)
- KpnI(495)
- BamHI(497)
- EcoRI(503)
- NruI(509)
- AsiSI(513)
- NotI(516)
- XhoI(524)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	StrepSUMOAvi

B42cl MCS region

XmnI BstBI SpeI KpnI BamHI EcoRI
 ...GCTCAGAAAATCGAATGGCACGAAGAAAAAGTTCTCGAACTAGTGGTACCGGATCCGAATTC
 Avi ...A Q K I E W H E E K V R T S G T G S E F

AsiSI NotI XhoI
 GCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 A I A A A L E H H H H...