



- XmnI(790)
- BstBI(794)
- KpnI(808)
- BamHI(810)
- EcoRI(816)
- NruI(822)
- AsiSI(826)
- NotI(829)
- XhoI(837)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	TrxHisSUMOstrep

B60cl MCS region

XmnI BstBI SpeI KpnI BamHI EcoRI
 ...CACCCGCAGTTTGAAAAAGGTGCGGAAAAAGTTCGAACTAGTGGTACCGGATCCGAATTC
 StrepII ...H P Q F E K G A E K V R T S G T G S E F

AsiSI NotI XhoI
 GCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 A I A A A L E H H H H...