



- XmnI(805)
- BstBI(809)
- KpnI(823)
- BamHI(825)
- EcoRI(831)
- NruI(837)
- AsiSI(841)
- NotI(844)
- XhoI(852)

Vector Features

Promoter	T7
Host Cell	<i>E. Coli</i>
Bacterial selection antibiotic	Ampicillin
Tag	Trx HisSumoAvi

B62cl MCS region

XmnI BstBI SpeI KpnI BamHI EcoRI
 ...GCTCAGAAAATCGAATGGCACGAAGAAAAAGTTTCGAACTAGTGGTACCGGATCCGAATTC
 Avi ...A Q K I E W H E E K V R T S G T G S E F

AsiSI NotI XhoI
 GCGATCGCGGCCGCACTCGAGCACCACCACCAC...
 A I A A A L E H H H H...