



- XmnI(661)
- BstBI(665)
- AgeI(669)
- KpnI(676)
- SmaI(678)
- EcoRI(681)
- XhoI(687)
- AsiSI(696)
- NotI(701)

Vector Features

Promoter	CMV
Host Cell	Mammalian
Bacterial selection antibiotic	Ampicillin
Mammalian selection marker	Yes
Tag	N-His

M01cl MCS region

...TACCATCACCATCACCATCACTCTAGATCAGAAGGAGTTCGAACCGGTACCCGGGAATTC
 ...Y H H H H H S R S E G V R T G T R E F

XhoI AsiSI NotI
 CTCGAGCGATCGCGGGCCGCAACCCAGCT...
 L E R S R G R