



BstBI(612)
AgeI(616)
KpnI(623)
SmaI(625)
EcoRI(628)
XhoI(634)
AsiSI(643)
NotI(648)

Vector Features

Promoter	CMV
Host Cell	Mammalian
Bacterial selection antibiotic	Ampicillin
Mammalian selection marker	Yes
Tag	N/A

M02cl MCS region

CMV...**CCGGGACCGATCCAGCCTCCG**GACTCTAGCCTAGATCTTCGAACCGGTACCCGGGAATTC

BstBI AgeI SmaI EcoRI

XhoI AsiSI NotI

CTCGAGCGATCGCGGGCCGCAACCCAGCT...