



- XmnI(718)
- BstBI(722)
- AgeI(726)
- KpnI(733)
- EcoRI(735)
- AsiSI(747)
- NotI(752)
- XhoI(760)
- PmeI(768)

### Vector Features

Promoter	CMV
Host Cell	Mammalian
Bacterial selection antibiotic	Ampicillin
Mammalian selection marker	Yes
Tag	Myc

### M43cl MCS region

Myc
XmnI
BstBI
AgeI
KpnI
EcoRI

...ATGGAGCAGAAACTCATCTCTGAAGAGGATCTGGAAAGAGTTCGAACCGGTACCGAATTC  
M E Q K L I S E E D L E R V R T G T E F

AsiSI
NotI
XhoI
PmeI

TCGCGATCGCGCGGCCGCGCTCGAGTTTAAACCTCGACT  
 S R S R G R A R V \*